



Bed Bug ALERT™ Bedbug Monitor Instructions for Use

UK Sole importer:

Bird-X UK & Ireland

Units 2c & 2d Littlehampton Marina
Ferry Road
Littlehampton
West Sussex BN17 5DS
United Kingdom
Tel: +44 (0)1903 715 631

Fax: +44 (0)1903 539 331
Email: sales@bird-x.co.uk
www.bird-x.co.uk and www.bird-x.ie

- Bed bugs hide in cracks, crevices and other small areas of safe harbourage such as between cushions and mattresses, making them hard to detect.
- This clear, safe to use, rigid plastic monitoring device slips between mattress and frame, or any place that you want to keep safe.
- The proprietary gel attracts and traps bed bugs if they are present.
- Each monitor will last up to three months if used in a clean environment.
- Use in: hotel/motel rooms, hospital wards, homes, dormitories, nursing homes, luggage racks/stores, cars, buses, lobbies, restaurant or anywhere with carpeting or soft furnishings..
- Results are possible within hours, though in most cases it will take several days.

EASY TO USE

- 1. Simply insert the business-card-size Bed Bug ALERT™ monitor between the mattress and box spring/divan, in the corners at the head or the bed, or into any crevice in any suspect area, i.e. under sofa cushions or behind furniture against skirting boards.
- 2. Check the Bed Bug ALERT™ daily and discard immediately in a sealed plastic bin bag should it contain a bed bug.
- 3. If ANY sign of bed bugs is found, contact your pest control service provider IMMEDIATELY to take pesticidal action to control and eliminate infestation immediately.

TO RE-ORDER

To re-order this product simply contact your pest control service provider or call Bird-X UK & Ireland on +44 (0)1903 715 631. Bed Bug ALERT monitors are sold in Packs of 12 monitors or cases of 144 monitors.