

MediK4n

Kills pathogens, coronaviruses,
including COVID-19.



MediK4n is a highly effective, powerful broad-spectrum sanitiser. MediK4n sanitises and helps prevent the spread of bacteria and viruses, it also eliminates unpleasant odours.

MediK4n allows you to quickly and easily sanitise and protect an entire space with just the push of a button.

A propellant inside the can sends out a fine aerosol mist of powerful, disinfecting active ingredients (Tested to BS EN 14476, 1276 and 1650) that disperse throughout the space. The aerosol gets to work on contact with every surface including difficult to reach areas, eliminating a wide range of germs, bacteria, yeasts and fungi. It is extremely effective against viruses and other harmful pathogens, proven to kill 99.999% of microorganisms.

DIRECTIONS MUST BE READ BEFORE USE:

- + Each can will treat an area of approximately 75m³ (L5m, W5m, H3m). Additional cans may be used to treat larger spaces. Lay the can(s) out as required for the area to be covered and ensure cans are evenly distributed
- + Ensure the area to be treated is vacated and no people, pets, aquatics, delicate plants or foliage are left in the vicinity. Where necessary display warning signs to prevent unauthorised access.
- + Prepare the area to be treated by:
 - + Ensuring all heat sources are removed / turned off and close windows
 - + For best performance, it is recommended that cleaning of surfaces is carried out prior to use
 - + Do not operate directly below smoke alarms
 - + Shake the can well for 20 seconds, remove lid, ensure the can is upright and press and lock the total release valve button until it clicks open
 - + When using multiple cans, work in sequence towards the exit route to avoid walking in the aerosol mist
 - + Exit the space and close the door
 - + Once activated it cannot be stopped until the can has fully dispensed (less than 5 minutes)
 - + Vacate the area for at least 1 hour or preferably overnight. Ventilate area for 15 mins (i.e. open windows/doors)

MediK4n

Kills pathogens, coronaviruses,
including COVID-19.



SALES POINTS:

- + Convenient - use where and when and as frequently as needed
- + Easy to use - Once activated Medik4n releases its contents within 2 minutes, creating an anti-viral fog which covers all contact surfaces, including those difficult to reach areas easily missed by traditional cleaning
- + Cost effective - no need for manual spraying or extended PPE
- + Flexible - single can for area up to 75/100m³ or multiple cans for larger areas
- + Safe to use - no residue and no wetting of surfaces
- + Powerful formula which remains active- Tested to BS 1276, BS 1650 & BS EN 14476
- + BRITISH MADE
- + Guaranteed prices - No increase for 12 month
- + Bespoke service
- + KILLS Coronaviruses including COVID-19